

**Title:** Long circulating tracer tailored for magnetic particle imaging

**Authors:**

Sitong Liu<sup>1,\*</sup>, Andreina Chiu-Lam<sup>1,\*</sup>, Angelie Rivera-Rodriguez<sup>2,\*</sup>, Ryan DeGroff<sup>2</sup>, Shehaab Savliwala<sup>1</sup>, Nicole Sarna<sup>2</sup>, Carlos Rinaldi<sup>1,2,+</sup>


**Affiliations:**

1. Department of Chemical Engineering, University of Florida, Gainesville, FL 32611
2. J. Crayton Pruitt Family Department of Biomedical Engineering, University of Florida, Gainesville, FL 32611-6131



\*These authors contributed equally to this work.

+Corresponding author



**Supplementary information**




**Figure S1. Custom 3D-printed sample bed for phantom and SPION sample MPI studies.** Dimensions of this custom-designed sample bed were derived from Magnetic Insight's 45mL sample tube design and modified to fit the custom sample holders that were designed and 3D-printed for this study.




**Figure S2. Vertical capillary tube holder designed and 3D-printed for LoD studies.** Custom sample holders were designed in Onshape and 3D printed using Formlabs Clear V4 resin with the Form 3 SLA 3D printer.




**Figure S3. Vertical 0.2 mL microcentrifuge tube holder designed and 3D-printed for MPI relax scan measurements.** The design accommodates 7 total tubes, with the entire lengths of 4 tubes able to be seen from a side view through exposed spaces. These 4 exposed spaces are 27.91 mm in vertical height with a center spacing of 22.4 mm.




**Figure S4. Rendering of two-part animal bed assembly designed and 3D-printed for MPI animal studies.** The top and bottom animal bed parts (left) assemble a whole animal bed (right) that can be installed on Magnetic Insight's MPI scanner. The top part was designed with a flat bottom for ease of transfer of the mice from MPI to CT. The bottom part houses an integrated tubing system for isoflurane administration during imaging.




**Figure S5. Magnetic characterization of ferucarbotran.** A) MH curve at 300 K for ferucarbotran in water. B) MH curves at 295, 305, 315K for ferucarbotran in TEGMA C) ZFC/FC at 10 Oe for ferucarbotran in TEGMA.




**Figure S6. Magnetic characterization of Synomag<sup>®</sup>-D.** A) MH curve at 300 K for Synomag<sup>®</sup>-D in water. B) MH curves at 295, 305, 315K for Synomag<sup>®</sup>-D in TEGMA C) ZFC/FC at 10 Oe for Synomag<sup>®</sup>-D in TEGMA.



**Figure S7. Physical, magnetic, and hydrodynamic size distribution histograms for A) ferucarbotran. B) Synomag®-D. C) RL-1A. D) RL-1B. E) RL-1C.**



**Figure S8. Dynamic magnetic susceptibility characterization of relaxation mechanism for all three tracers. A) ferucarbotran. B) Synomag<sup>®</sup>-D. C) RL-1C. in water (red circle) or water and glycerol mixture (blue triangle, 65 wt% of glycerol). Dash line indicating peak frequency position for each particle assuming Brownian Relaxation dominant, based on their hydrodynamic diameter.**



**Figure S9. PSF for several RL-1 tracer batches demonstrate reproducibility in MPI performance.**




Figure S10. Representative 2D MPI z-channel dilution series for commercially available tracers and RL-1C. MPI scans were acquired in high-sensitivity.





Figure S11. Mean signal to noise ratio (mSNR) as a function of tracer mass.



**Figure S12. Collage showing representative MPI images to each time point for one mouse for RL-1C tracer.** The coregistered optical image with MPI is for reference to the signals observed in the MPI scans at different times. The ROIs for the heart and liver/spleen were all the same size across all times.

**Movie S1. Video of mouse injected with RL-1C tracer 1 hour after administration.** The video corresponds to coregistered CT image with 3D MPI scan.

**Movie S2. Video of mouse injected with RL-1C tracer 24 hours after administration.** The video corresponds to coregistered CT image with 3D MPI scan.

**Movie S3. Video of mouse injected with ferucarbotran tracer 24 hours after administration.** The video corresponds to coregistered CT image with 3D MPI scan.

**Movie S4. Video of mouse injected with Synomag<sup>®</sup>-D tracer 24 hours after administration.** The video corresponds to coregistered CT image with 3D MPI scan.